




The Homogenizer Company™



Bead Ruptor 4

MINI BEAD MILL HOMOGENIZER

Built upon the innovative design and unsurpassed quality of Omni's Bead Ruptor family.

***Powerful performance
in a compact design***

Bead Ruptor 4

MINI BEAD MILL HOMOGENIZER

ORDER NO. 25-010



Processing Range:

25 μ L to 6 mL

Power Rating:

100-230 VAC/50-60 Hz

Dimensions:

8.5" W x 11.5" D x 10" H

Weight: 16.5 lbs. (7.48 kg)

Run Time: 1 - 5 mins.

Programs: 1

Cycles: 1

Performance Range:

30 g - 150 g

Efficiently lyse:

- > Human tissues
- > Animal tissues
- > Plant & food
- > Micro-organisms

Designed for personal use

- > Broad performance range — gentle mixing to high force homogenization
- > Small footprint — requires minimal bench space
- > Holds 4 x 0.5 mL / 1.5 mL / 2 mL or 1 x 7 mL tubes

Patent Pending Snap In Tube Carriage



Omni Bead Ruptor 4 is ideally suited for:

- > Homogenizing up to 4 samples in seconds
- > Reproducible extraction of proteins and nucleic acids
- > Simultaneously processing 4 x 0.5 mL, 4 x 1.5 mL, 4 x 2 mL or 1 x 7 mL samples



For best results, use with Omni pre-filled bead kits



Omni International, Inc. www.omni-inc.com
935C Cobb Place Blvd. Kennesaw, GA 30144
800.776.4431 770.421.0058

