




The Homogenizer Company™



Bead Ruptor 24

BEAD MILL HOMOGENIZER

*Quickly disrupt
the toughest samples*

**Omni Bead Ruptor 24
is ideally suited for:**

- > DNA/RNA extraction
- > Tissue homogenization
- > Enzyme isolation
- > Protein purification
- > Food safety testing
- > Drug detections
- > qPCR, RT-PCR, electrophoresis,
mass spectrometry

Bead Ruptor 24

Bead Mill Accessories



96 Strip Well Tube Carriage Kit

ORDER NO. 19-010-330

- > 96 x 1.2 mL strip well tubes
- > Includes tube carriage that holds the 1.2 mL strip well tubes firmly in place during processing.



1.5 mL Microtube Carriage Kit

ORDER NO. 19-010-310S

- > 24 x 1.5 mL microtubes
- > Includes Bead Ruptor Finger Plate that holds the 1.5 mL microtubes firmly in place during processing.



2 mL Tube Carriage Kit

ORDER NO. 19-010-310

- > 24 x 2 mL or 24 x 0.5 mL skirted tubes
- > Includes Bead Ruptor Finger Plate that holds the 2 mL tubes firmly in place during processing.



7 mL Tube Carriage Kit

ORDER NO. 19-010-320

- > 12 x 7 mL tubes
- > Includes Bead Ruptor Finger Plate that holds the 7 mL tubes firmly in place during processing.



30 mL Tube Carriage Kit

ORDER NO. 19-345-035

- > 6 x 30 mL tubes
- > Patent pending snap-in tube carriage for 30 mL samples. Kit includes tube carriage and lock ring.



50 mL Tube Carriage Kit

ORDER NO. 19-345-050

- > 3 x 50 mL tubes
- > Patent pending snap-in tube carriage for 50 mL samples. Kit includes tube carriage and lock ring.

ORDER NO. 19-040

> Specifications

Processing Range: 25 μ L to 50 mL

Power Rating: 100-230 VAC/
50-60 Hz

Dimensions: 15" W x 17" D x 13" H

Weight: 29 kg (65 lbs.)

Compliance: CE

Run Time: 1 s - 9:59 m

Programs: 99

Cycles: 1-10

Performance Range:
0.8 m/s - 8 m/s
in increments of 0.15 m/s



Omni International, Inc. www.omni-inc.com
935C Cobb Place Blvd. Kennesaw, GA 30144
800.776.4431 770.421.0058