

## COMBO OFFER PRICE ₹ 92,500

### ENDURO™ Vertical Gel With Power Supply

#### ENDURO™ VE10 VERTICAL GEL SYSTEM

Labnet's ENDURO™ Modular Mini Vertical Gel System for protein electrophoresis is ideal for screening new samples and evaluating sample preparation conditions.

##### SPECIFICATIONS

Gel (W x L)	10 x 10 cm
Gel Capacity	4 gels
Unit (W x D x H)	19 x 13 x 15 cm
Buffer Volume	250 mL
Max. Samples	20 per gel



#### ENDURO™ POWER SUPPLY 300V



The ENDURO™ 300V Power Supply is ideal for general laboratory use especially for horizontal gel electrophoresis and mini-vertical PAGE gels. This high value power supply provides operating choices that include constant voltage or constant amperage with automatic crossover, no load detection, sudden load change and over voltage protection. It comes with a timer and

allows for timed or continuous operation. All operating parameters are set using the membrane keypad and displayed digitally on the backlit LCD. Safety features are incorporated into the power supply systems. The units can also be programmed to restart in the event of a power disruption. This comes with to run up to 4 gel tanks simultaneously.

##### SPECIFICATIONS

	Model 300V
Output voltage range / increments	2-300V / 1V
Output current range / increments	4-500 mA / 1 mA
Output power range / increments	90 watts max. / na
Operating constant modes	Voltage or current
Timer	1-999 minutes or continuous
Programmable	No
Dimensions (W x D x H)	19 x 25 x 8 cm
Weight	2.5 kg
Input Voltage	120-230V

##### ORDERING INFORMATION

Catalog No.	Description
E2010-PA	ENDURO PAGE system, includes PAGE insert, buffer tank with leads and cooling pack
E0303	ENDURO 300 volt Power Supply

**Offer Ends on  
15<sup>th</sup> Oct. 2016**

• *Taxes extra as applicable*

• *Contact us : [sales@gingerscience.co.in](mailto:sales@gingerscience.co.in) | [yogendra@gingerscience.co.in](mailto:yogendra@gingerscience.co.in)*

Order to be placed to :

**Ginger Science Pvt. Ltd.**

LGF-6-8, Jinkushal Suri Complex

230-A, Masjid Moth, South Ex-II | New Delhi – 110049

Cell No: +91-9650220266 | Tel No : 011-26263236 /26253236

[www.gingerscience.co.in](http://www.gingerscience.co.in)