

Microcentrifuge Tubes

Corning offers two styles of microcentrifuge tubes: traditional snap cap tubes for quick access or screw cap tubes for greater sealing security.

Snap Cap Polypropylene Microcentrifuge Tubes



Microcentrifuge Tubes

- ▶ Costar® microcentrifuge tubes are certified RNase-/DNase-free
- ▶ Supplied nonsterile and are autoclavable
- ▶ External graduations and frosted writing spot for easy sample identification
- ▶ Positive seal design allows for repeated opening and closing
- ▶ Flat cap surface for convenient labeling
- ▶ Withstands a maximum RCF of 17,000 xg
- ▶ Costar low binding microcentrifuge tubes feature a bonded polymer technology that reduces protein and nucleic acid binding, resulting in better sample recovery

Snap Cap Polypropylene Microcentrifuge Tubes Ordering Information

Cat. No.	Volume (mL)	Color	Qty/Pk	Qty/Cs
Snap Cap Microcentrifuge Tubes				
3208	0.65	Natural	500	1,000
3209*	0.65	Rainbow*	200	1,000
3620	1.7	Natural	500	500
3621	1.7	Natural	500	5,000
3622*	1.7	Rainbow*	100	500
3213	2.0	Natural	500	1,000
Low Binding Snap Cap Microcentrifuge Tubes Ordering Information				
3206	0.65	Natural	500	500
3207	1.7	Natural	250	250


*Rainbow pack includes one bag each of blue, green, yellow, red, and purple tubes.

Screw Cap Polypropylene Microcentrifuge Tubes




Microcentrifuge Tubes

- ▶ Corning® polypropylene microcentrifuge tubes feature screw caps that provide a tight secure seal
- ▶ Choice of attached cap with silicone O-ring or unattached rim seal cap
- ▶ Withstands a maximum 20,000 RCF
- ▶ Sterile



▶ Attached loop cap allows for optimum one-handed convenience. Silicone O-ring gasket provides a snug seal, safeguarding samples against leakage.



▶ Easy-to-use unattached rim seal cap design twists on or off in a single turn.

Screw Cap Polypropylene Microcentrifuge Tubes Ordering Information

Cat. No.	Volume (mL)	Cap Style	O-ring	Self Standing	Qty/Cs
430909	1.5	Attached	Yes	No	500
430915	2.0	Attached	Yes	Yes	500
430917	2.0	Unattached	No	Yes	500

Can't Find the Products You Need?

Call us at 800.492.1110 or 978.442.2200 to find out about our new range of Axygen® brand Microtubes.