



**Buy 1 and Get 1 FREE!**

OFFER PRICE – ₹ 15,800 for Lambda Plus Single Channel Pipettor  
 ₹ 54,000 for Lambda Plus 8 Channel Pipettor  
 ₹ 63,000 for Lambda Plus 12 Channel Pipettor



Lambda Plus Single-Channel Pipettor

### Corning Lambda Plus Single-Channel, 8-Channel and 12-Channel Pipettors

Lambda Plus Pipettors have been engineered to provide the highest level of accuracy, reproducibility, and comfort. Available in single-, eight- and twelve-channel configurations, all Lambda Plus Pipettors feature smooth plunger movement and extremely low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button, and a locking mechanism prevents accidental change.

The Lambda Plus Pipettors have been designed to maintain their calibration, even after repeated use and autoclaving. A key is included with each pipettor for easy in-house calibration.

- Factory calibrated to strict ISO 8655 standards
- 3-year warranty
- Fully autoclavable and UV resistant
- Designed for universal fit with all common brands of pipette tips
- Adjustable tip ejector



Lambda Plus 8-Channel Pipettor

### Lambda Plus Single-Channel Pipettor Ordering Information

Cat. No	Volume Range (µL)	Accuracy	Precision	Tip (µL)	Qty/Cs
4070	0.2-2.0	±1.2 to ±1%	<6.0 to <0%	10	1
4071	0.5-10	±4.0 to ±0%	<4.0 to <0%	10	1
4072	2-20	±4.0 to ±0%	<3.0 to <0%	200	1
4073	10-100	±1.6 to ±0%	<0.8 to <0%	200	1
4074	20-200	±1.2 to ±0%	<0.6 to <0%	200	1
4075	100-1,000	±1.0 to ±0.6%	<0.45 to <0.2%	1,000	1



Lambda Plus 12-Channel Pipettor

### Lambda™ Plus 8-Channel and 12-Channel Pipettor Ordering Information

Cat. No.	Volume Range (µL)	Channels	Accuracy	Precision	Tip (µL)	Qty/Cs
4080	1-10	8	±8.0 to ±2.0%	<6.0 to <1.2%	10	1
4081	5-50	8	±4.0 to ±1.6%	<2.5 to <0.6%	200	1
4082	20-200	8	±3.0 to ±1.0%	<1.5 to <0.6%	200	1
4083	50-300	8	±1.6 to ±1.0%	<1.5 to <0.6%	200/300	1
4084	1-10	12	±8.0 to ±2.0%	<6.0 to <1.2%	10	1
4085	5-50	12	±4.0 to ±1.6%	<2.5 to <0.6%	200	1
4086	20-200	12	±3.0 to ±1.0%	<1.5 to <0.6%	200	1
4087	50-300	12	±1.6 to ±1.0%	<1.5 to <0.6%	200/300	1

- Taxes extra as applicable
- Contact us : [sales@gingerscience.co.in](mailto:sales@gingerscience.co.in) | [yogendra@gingerscience.co.in](mailto:yogendra@gingerscience.co.in)

Order to be placed to :

**Ginger Science Pvt. Ltd.**

LGF-6-8, Jinkushal Suri Complex

230-A, Masjid Moth, South Ex-II | New Delhi – 110049

Cell No: +91-9650220266 | Tel No : 011-26263236 / 26253236

**Offer Ends on 31<sup>st</sup> Oct. 2016**