



8- and 12-Pette Multichannel Pipettors

8-Pette® and 12-Pette Multichannel Pipettors

These pipettors feature a unique, ergonomic trigger-style aspiration and dispense control mechanism designed to reduce thumb fatigue during repetitive pipetting. Volume range is 20 to 200 µL and can be adjusted with a vernier-scale spindle. Autoclavable.

8-Pette and 12-Pette Multichannel Pipettors Ordering Information

Cat. No.	Volume Range (µL)	Channels	Qty/Cs
4880	20-200	12	1
4888	20-200	8	1



Stripettor Plus Pipetting Controller with pipet

Corning® Stripettor™ Plus Pipetting Controller

The new Corning Stripettor Plus is comfortable, lightweight, and designed for use with most brands of serological pipets. Speed range is adjustable via a control wheel on the back of the instrument, and aspirating and dispensing can be finely controlled by the pressure-sensitive buttons. A convenient wall hanger, charging stand, and spare filters are included with the Stripettor Plus, and when fully charged, the high capacity Lithium-Ion battery allows continuous pipetting for up to 4.5 hours. A battery indicator light lets the user know when it's time to recharge; the Stripettor Plus can be used while charging. The charging stand allows storage on its back, even with pipet inserted.

- ▶ Ergonomic, comfortable design
- ▶ Wall hanger, charging stand, charging power supply, and spare filters included
- ▶ Variable speed range and pressure sensitive control buttons
- ▶ Settings for blow out or slow “gravity” dispensing
- ▶ Compatible with pipets from 1 mL to 100 mL
- ▶ Includes hydrophobic, PTFE 0.2 µm filters
- ▶ Fast aspiration and dispensing speed at highest setting: 25 mL in 2.5 seconds
- ▶ Li-Ion battery provides 4.5 hours of continuous use



Stripettor Plus Pipetting Controller in charging position

Stripettor Plus Pipetting Controller and Accessories Ordering Information

Cat. No.	Description	Voltage	Qty/Cs
4090	Stripettor Plus with three 0.2 µm filters, wall hanger and charger	120V, 50/60Hz	1
4091	Stripettor Plus with three 0.2 µm filters, wall hanger and charger	230V, 50/60Hz, EU Plug	1
4092	Stripettor Plus with three 0.2 µm filters, wall hanger and charger	230V, 50/60Hz, UK Plug	1
4093	Stripettor Plus with three 0.2 µm filters, wall hanger and charger	100V, 50/60Hz, JP Plug	1
4097	Stripettor Plus with three 0.2 µm filters, wall hanger, without charger	230V, 50/60Hz, KR Plug	1
4098	Stripettor Plus with three 0.2 µm filters, wall hanger and charger	230V, 50/60Hz, AU Plug	1

EU = Europe, UK = United Kingdom, JP = Japan, KR = Korea, AU = Australia

Accessories

4094	PTFE autoclavable filter, 0.2 µm	5
4095	PTFE autoclavable filter, 0.45 µm	5
4096	Replacement silicone pipet holder	1



Stripettor Plus Pipetting Controller in storage position

Corning Aspirator Device and Accessories

Aspirator and Accessories Ordering Information

Cat. No.	Description
4930	Aspirator device (includes hand piece, grommet for accessory attachment, and single-channel adapter for use with disposable pipet tips)
4931	8-channel adapter for use with disposable pipet tips
9016	Costar® aspirating pipet, without graduations, unplugged, sterile, paper/plastic wrapped



Aspirator Device