

ProBlot™ Rocker 25 and 25XL

The Labnet compact rockers, *ProBlot 25* and *25XL*, are ideal for staining, hybridization and general mixing. Their basic design and economical price makes them a great value for any laboratory.

Rocking speed is adjustable to provide the gentle motion required for staining fragile gels as well as the more vigorous action necessary for washing blots. The rocking movement is fixed at 7°.

The platform surfaces are covered with non-slip rubber pads to prevent sample containers from moving during operation. Elastic tie downs are available to secure bags, boxes, tubes and other containers in place. The *ProBlot 25* platform measures 26 x 20 cm, while the *ProBlot 25XL* platform is significantly larger at 30 x 30 cm. For increasing the work area without sacrificing valuable bench space, double platform models can be purchased.

The *ProBlot 25* and *25XL* Rockers can be used safely in temperature controlled environments up to 70°C. Their low profile and compact size make them perfect for use in incubators.

- Designed for staining, hybridization, washing, etc.
- Variable rocking speed
- Compact size, low profile
- Large, corrosion resistant stainless steel platforms
- Safe for coldroom or incubator use



CE

SPECIFICATIONS

Speed settings	5 to 120 tilts/minute
Motion	Rocking
Tilt range	Fixed, 7°
Clearance between platforms	3.5 in/9 cm
Maximum load	11 lb/5 kg
Ambient operating range	+4° to 70°C
Dimensions (W x D x H)	
ProBlot 25	11.4 x 12.6 x 5.7 in/ 29 x 32 x 14.5 cm
ProBlot 25XL	12.8 x 14.2 x 5.7 in/32.5 x 36 x 14.5 cm
Weight	11.7 lb/5.3 kg
Electrical	230V~, 50/60 Hz 120V~, 50/60 Hz

ORDERING INFORMATION

Catalog No.	Description
S2025-230V*	<i>ProBlot Rocker 25</i> with single platform (26 x 20 cm), 230V
S2025-D-230V*	<i>ProBlot Rocker 25D</i> with double platform (26 x 20 cm), 230V
S2025-XL-230V*	<i>ProBlot Rocker 25XL</i> with large single platform (30 x 30 cm), 230V
S2025-XLD-230V*	<i>ProBlot Rocker 25XLD</i> with large double platform (30 x 30 cm), 230V
S2025-T	Elastic tie downs, pack of 4
S2025-ST	Gel staining tray, pack of 4

*230V product includes both EU and UK power cords. For 120V units, remove -230V from end of catalog number.