

EQUIPMENT



Low Speed Orbital Shaker

Corning® LSE™ Low Speed Orbital Shaker

The Corning LSE Low Speed Orbital Shaker includes a 30 x 30 cm platform and is ideal for staining and destaining gels, washing blots and general mixing applications. Speed can be set from 3 to 60 rpm and the timer can be set for up to 2 hours or for continuous operation.

- ▶ Ideal for staining and washing of fragile gels
- ▶ Gentle motion prevents foaming of liquids
- ▶ Simple analog control with easy-turn knobs
- ▶ 19 mm orbit for blotting, gel staining and general mixing
- ▶ Safe for cold room and incubator use

Corning LSE Low Speed Orbital Shaker Ordering Information

	Cat. No.	Description	Voltage	Plug Type	Speed (rpm)	Weight (kg)	Dimensions (cm)	Qty/Cs
New	6780-FP	LSE Low Speed Orbital Shaker	120V, 50/60Hz	US	3-60	5	25.5 x 32 x 16	1
New	6781-FP	LSE Low Speed Orbital Shaker	230V, 50/60Hz	EU	3-60	5	25.5 x 32 x 16	1
New	6782-FP	LSE Low Speed Orbital Shaker	230V, 50/60Hz	UK	3-60	5	25.5 x 32 x 16	1



Orbital Shaker with options clamps and flasks

Corning LSE Orbital Shaker

The Corning LSE Orbital Shaker meets a variety of shaking and mixing applications. A powerful brushless motor and counter-balanced shaker mechanism provide smooth operation for years of trouble-free use. Different platforms (order separately) are available and are easily interchangeable.

- ▶ Two stainless steel platform options
- ▶ 19 mm orbit for aeration of cultures and staining gels
- ▶ Large digital display with single-knob control
- ▶ Digital control of shaking speed and time
- ▶ Safe for cold room or incubator use
- ▶ Maximum capacity: 4 x 1L, 5 x 500 mL, 9 x 250 mL or 16 x 125 mL

Corning LSE Orbital Shaker Ordering Information

	Cat. No.	Description	Voltage	Plug Type	Speed (rpm)	Weight (kg)	Dimensions (cm)	Qty/Cs
New	6780-NP	LSE Orbital Shaker	120V, 50/60Hz	US	20-300	7.3	25.6 x 31.8 x 12.7	1
New	6781-NP	LSE Orbital Shaker	230V, 50/60Hz	EU	20-300	7.3	25.6 x 31.8 x 12.7	1
New	6782-NP	LSE Orbital Shaker	230V, 50/60Hz	UK	20-300	7.3	25.6 x 31.8 x 12.7	1

Accessories

New	480108	Flat Platform with Rubber Mat						1
New	480109	Flask Clamp Platform (without clamps)						1
New	480110	Clamp for 25 mL flask						1
New	480111	Clamp for 50 mL flask						1
New	480112	Clamp for 125 mL flask						1
New	480113	Clamp for 250 mL flask						1
New	480114	Clamp for 500 mL flask						1
New	480115	Clamp for 1L flask						1
New	431405	Erlenmeyer flask, baffled, 125 mL, vent cap, sterile						50
New	431407	Erlenmeyer flask, baffled, 250 mL, vent cap, sterile						50
New	431401	Erlenmeyer flask, baffled, 500 mL, vent cap, sterile						25
New	431403	Erlenmeyer flask, baffled, 1L, vent cap, sterile						25



Polycarbonate Erlenmeyer Flasks, Baffled

For a full listing of Corning flasks, please visit www.corning.com/lifesciences.